# BROEN-LAB EMERGENCY SHOWER SYSTEMS

## General instructions about installation, maintenance and operation

94G0053 - 14.01.2025

These instructions are based upon the most common international norms, EN15154 1, 2 & 5 and ANSI Z 358.1. Remember to follow the norm, that applies with your geographical area.

EN15154 1, 2 & 5 is set by the European committee for Standardization. Part 1 covers "Plumbed-in body showers for laboratories", Part 2 covers "Plumbed-in eye wash units" and Part 5 covers "Water overhead body showers for sites other than laboratories".

ANSI Z 358.1 is set by the American National Standard Institute. It covers Emergency Eyewash and Shower Equipment.

EN15154-1, 2 & 5	English	Page 2
EN15154-1, 2 & 5	Danish	Page 4
EN15154-1, 2 & 5	German	Page 6
EN15154-1, 2 & 5	French	Page 8
EN15154-1, 2 & 5	Spanish	Page 10
ANSI Z 358.1	English	Page 12

## EN15154-1, 2 & 5 INSTALLATION OF EMERGENCY SHOWERS

BROEN-LAB Emergency showers are used as first aid in accidents causing burns, scalds, chemical contamination etc. on the human body or eye. The wording "Emergency shower" refers in the following to the complete range of body, eye and hand-held showers.

**Placement:** Less than 20 meter without obstacles from potential injury or max. 10 seconds from potential injury.

**Open space:** There must be an open space under the body shower of at least 400 mm in radius. See picture 1.

A marking on the floor is advisable to remind the workers to keep the area clear.

**Water quality:** Mains water only. Drinking water or equivalent quality.

**Minimum dynamic pressue:** As indicated on sticker on product. Must be obtained during operation.

**Flow rate:** Body shower 60 l/min  
Eye shower 6 l/min per shower head

**Water available for:** Min. 15 minutes

**Water temperature:** Recommended between 15 - 37° C (idealy 20 - 25° C)

**Spray pattern:** At a distance of 700 mm below the shower head, 50 (± 10) % of the volume of water delivered shall fall in a circle with a radius of 200 mm. See picture 2.

**Safety signs:** Signs in accordance to ISO 3864-1 is included in the package. See our recommendations for mounting height in picture 3.


For Australia signs according to AS 4775 / AS 1319 are included as well.

**Mounting height:** See table 1 and picture 3.


1. Eye shower, distance from floor	1000 ± 200 mm
2. Distance between sign and product	Max. 250 mm *
3. Sign on the wall, distance from floor	Max. 2000 mm *
4. First aid plan, distance from floor	Max. 1500 mm *
5. Triangular sign, distance from product	Max. 250 mm *
6. Eye wash station, distance from floor	Max. 1700 mm *
7. Service sign, distance from floor	Max. 1400 mm *
8. Pull rod, distance from floor	Max. 1750 mm
9. Body Shower, distance from floor	2200 ± 100 mm
10. Sign placed	Center of pullrod

Table 1: Mounting height


\* Recommendation from BROEN-LAB

Picture 1: Open space

Picture 2: Spray pattern

Picture 3: Mounting height

## TEST AND MAINTENANCE OF EMERGENCY SHOWERS

According to the EN15154-1 the manufacturer shall supply a method and frequency for testing. BROEN-LAB recommends following the EN15154-5, that says that emergency showers shall be tested at least monthly or more according to national regulations if required.

**Flow rate:** Test if the flow rate is according to the description on the left; 60 l/min for body showers and 6 l/min per shower head for eye shower.  
A BROEN-LAB test trolley or measuring cylinder can be used for this purpose.

**Water quality:** Make a visual test of the water to see, if the water is clear, and therefore in accordance with the requirement for drinking water quality or equal.

**Spray pattern:** The water must be distributed in a conical shape in accordance with picture 2.  
A BROEN-LAB test trolley can be used for this purpose.

To comply with EN15154-5 the test needs to be documented.

**Maintenance:** BROEN-LAB recommends emptying the hand-held eye shower head as much as possible after using the eye shower. To do this, the handheld eye shower is held upside down when it is closed.

To avoid limescale deposits, the outer wetted parts (shower head, rubber cap and handle if necessary) should be wiped dry with a clean cloth. Limescale deposits on the shower head and other parts can thus be greatly reduced. This can extend the service life of the shower head until it is decalcified or replaced.

BROEN-LAB recommends to only use cleaning material which is explicitly provided for this type of application, preferably a mild, citric acid-based cleaner. We do not recommend mixing different cleaning agents. After cleaning, rinse thoroughly with clean water to remove any cleaner residue.

## OPERATION OF EMERGENCY SHOWERS

**Activation:** Water will flow automatically once the valve is opened.

Body shower	Pull down the pull rod to open the valve fully
Eye wash	Push activation arm back to open the valve fully
Hand-held eye shower	Squeeze the red handle to open the valve fully

**Plan of action:** If an accident arise, follow the steps below

- Start flushing immediately - every second counts
- Remove contaminated clothing, which is not stuck to skin
- In case of an eye injury, keep the eyelid open when continuing flushing
- Call for medical assistance and follow their instructions

**Inform employees:** All employees should have knowledge about

- How to prevent accidents/injuries
- How to react, should an accident arise
- The location of the emergency showers
- How to operate the emergency showers
- How to help injured person in the best way

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## INSTALLATION OF EMERGENCY SHOWERS

BROEN-LAB emergency showers are used as first aid in accidents causing burns, scalds, chemical contamination etc. on the human body or eye. The wording "Emergency shower" refers in the following to the complete range of body, eye and hand-held eye showers.

**Placement:** Approx 55 ft (16.8 m) to be reached within 10 seconds when walking

**Open space:** There must be an open space under the body shower of at least 40.6 cm (16 in.) in radius. See picture 1.

A marking on the floor is advisable to remind the worker to keep the area clear.

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**Minimum dynamic pressure:** As indicated on sticker on product. Must be obtained during operation.

**Flow rate:**

Body shower	75.7 l/min (20.0 gpm)
Eye wash	1.5 l/min (0.4 gpm)
Eye/face wash	11,4 l/min (3.0 gpm) - if it has two heads


**Water available for:** Min. 15 minutes

**Water temperature:** Recommended between 16 - 38°C (60 - 100°F)


**Spray pattern:** At a distance of 152.4 cm (60 in.) above the floor, the water shall fall in a circle with a diameter of 50.8 cm (20 in.) See picture 2.

**Safety signs:** Signs in accordance to ISO 3864-1 is included in the package. See our recommendations mounting heights in picture 3.

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Picture 1: Open space


Picture 2: Spray pattern

1. Eye shower, distance from floor	83.8 - 134.6 cm (33 - 53 in)
2. Distance between sign and product	Max. 25.4 cm (10 in) *
3. Sign on the wall, distance from floor	Max. 200,6 cm (79 in) *
4. First aid plan, distance from floor	Max. 149.9 cm (59 in) *
5. Triangular sign, distance from product	Max. 25.4 cm (10 in) *
6. Eye wash station, distance from floor	Max. 175.3 cm (69 in) *
7. Service sign, distance from floor	Max. 139.7 cm (55 in) *
8. Pull rod, distance from floor	Max. 175.3 cm (69 in) *
9. Body Shower, distance from floor	208.3 - 243.8 cm (82 - 96 in)
10. Sign placed	Center of pullrod

Table 1: Mounting height

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Emergency safety showers should be tested at least every week or more according to national regulations if required.

- Flow rate:** Using a flowmeter or other means, determine that the rate of flow according to the recommendations; 75.7 l/min (20 gpm) for body shower, and 1.5 l/m (0.4 gpm) for eye wash. A BROEN-LAB test trolley or measuring cylinder can be used for this purpose.
- Water quality:** Make a visual test of the water to determine, if the water is clear, and therefore in accordance with the requirement for drinking water quality or equal.
- Temperature:** Using a temperature gauge or other means, determine that the flushing fluid is tepid.
- Spray pattern:** Measure the diameter of the spray pattern in accordance with picture 2. A BROEN-LAB test trolley can be used for this purpose.
- Maintenance:** BROEN-LAB recommends emptying the hand-held eye shower head as much as possible after using the eye shower. To do this, the hand-held eye shower is held upside down when it is closed.
- To avoid limescale deposits, the outer wetted parts (shower head, rubber cap and handle if necessary) should be wiped dry with a clean cloth. Limescale deposits on the shower head and other parts can thus be greatly reduced. This can extend the service life of the shower head until it is decalcified or replaced.
- BROEN-LAB recommends to only use cleaning material which is explicitly provided for this type of application, preferably a mild, citric acid-based cleaner. We do not recommend mixing different cleaning agents. After cleaning, rinse thoroughly with clean water to remove any cleaner residue.

## OPERATION OF EMERGENCY SHOWERS

- Activation:** Water will flow automatically once the valve is opened.
- |                      |  |
|----------------------|--|
| Body shower          | Pull down the pull rod to open the valve fully   |
| Eye wash             | Push activation arm back to open the valve fully |
| Hand-held eye shower | Squeeze the red handle to open the valve fully   |

- Plan of action:** If an accident arise, follow the steps below
- Start flushing immediately - every second counts
  - Remove contaminated clothing, which is not stuck to skin
  - In case of an eye injury, keep the eyelid open when continuing flushing
  - Call for medical assistance and follow their instructions

- Inform employees:** All employees should have knowledge about
- How to prevent accidents/injuries
  - How to react, should an accident arise
  - The location of the emergency showers
  - How to operate the emergency showers
  - How to help injured person in the best way

W

A wellspring  
of certainty

**BROEN-LAB**

### **Warning!**

Before the installation of your BROEN-LAB emergency shower, all debris which may cause leaks must be thoroughly flushed from the system. Should a fault appear in this product, and the system was found not to have been flushed out, the warranty will be void.

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### **Contact a specialist**

– and find out more about what  
BROEN-LAB can do for your business

## **BROEN-LAB**

### **BROEN-LAB A/S**

Drejervaenget 2  
5610 Assens  
Denmark

Tel. +45 6376 6376

Fax. +45 6471 2476

lab@broen-lab.com

broen-lab.com