

Yearly safety training

In general the safety training takes place at your company and in immediate extension of one of our service visits.

You are able to book one training each year.

A safety training with BROEN-LAB includes:

Introduction to EN15154 and other relevant safety standards that we test in accordance to

Demonstration of how we test the emergency showers and what we test for

Installation of emergency showers, e.g. installation height relevant to accessibility

Identification of the locations of the emergency breakers and signage in relation to them

Description of the hazards present on the site (e.g. steam, chemicals, heat, dust, fire) and how they must be handled

Description of which hazards are near the various emergency shower stations

Accessibility to the shower – including free space under the shower and between the shower and potential danger

Activation of body showers

Activation of eyewashes

First aid in connection with an accident – e.g. flushing times and use of pH-neutralizing wash bottles

Colleague help and first aid in relation to flushing and partial flushing

If you have specific need/request it can be included in the drawing